

REPORT DOCUMENTATION PAGE

AFRL-SR-AR-TR-04-

Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering the required data, completing and reviewing this collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports (0704-0188), 4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to any penalty for failing to comply with a collection of information if it does not have a valid OMB control number. PLEASE DO NOT RETURN YOUR FORM TO THE ABOVE ADDRESS.

1. REPORT DATE (DD-MM-YYYY) Dec 10 th 2003		2. REPORT TYPE Final		3. DATES COVERED (From - To) Apr 01 2002 to Apr 01 2003	
4. TITLE AND SUBTITLE Synthesis of Oxide Fibers Using a Fiber Drawing Machine				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER F49620-02-1-0285	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S) Professor Waltraud M. Kriven 20040130 081				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Materials Science and Engineering Dept., University of Illinois at Urbana-Champaign, 1304 W. Green St., Urbana IL 61801 Tel: (217) 333 5258 Fax: (217) 333 2736 Email: kriven@uiuc.edu				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING / MONITORING AGENCY NAME(S) AND ADDRESS(ES) Dr. Joan Fuller Ceramic and Non-Metallic Materials Program Manager Directorate of Aerospace and Materials Engineering Air Force Office of Scientific Research 4015 Wilson Blvd, room 713 Arlington, VA 22203-1977 Ph. (703) 696 7236				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION / AVAILABILITY STATEMENT Unlimited					
13. SUPPLEMENTARY NOTES					
14. ABSTRACT The Marksman Fiber Drawing Machine, (Chemat Technology Inc., Northridge, CA) is a research scale, positive displacement, piston-type, fiber extruder. The machine was custom made to extrude various compositions of fibers from polymeric precursors specifically designed with sufficient flexibility to cover anticipated research requirements. It has a capacity of 10-150 grams, offers operational temperatures up to 500°C, a 10 – 10,000 m/min take up capacity, and 2400 Watts heater capacity. It allows the continuous drawing of fiber while the diameter of the spinnerette hole can be designed to fit application requirements. Standard spinnerette orifice diameters of 200 nm, 150 nm and 100 nm were supplied. In addition, an orifice diameter of 30 nm was custom made for research use. The machine for UIUC was custom made with an integrated high temperature furnace to allow the fiber to be dried and sintered at elevated temperatures, in three zones, up to 1,200°C with programmable controls.					
15. SUBJECT TERMS Fiber extruder, polycrystalline powders, 200nm to 30 nm orifice					
16. SECURITY CLASSIFICATION OF: none			17. LIMITATION OF ABSTRACT	18. NUMBER OF PAGES 4	19a. NAME OF RESPONSIBLE PERSON Prof. Waltraud Kriven
a. REPORT	b. ABSTRACT	c. THIS PAGE			19b. TELEPHONE NUMBER (include area code) (217) 333 5258

Description of Instrument

The Marksman Fiber Drawing Machine, (Chemat Technology Inc., Northridge, CA) is a research scale, positive displacement, piston-type, fiber extruder. The machine was custom made to extrude various compositions of fibers from polymeric precursors specifically designed with sufficient flexibility to cover anticipated research requirements. It has a capacity of 10-150 grams, offers operational temperatures up to 500°C, a 10 – 10,000 m/min take up capacity, and 2400 Watts heater capacity. It allows the continuous drawing of fiber while the diameter of the spinnerette hole can be designed to fit application requirements. Standard spinnerette orifice diameters of 200 nm, 150 nm and 100 nm were supplied. In addition, an orifice diameter of 30 nm was custom made for our research use. The machine for UIUC was custom made with an integrated high temperature furnace to allow the fiber to be dried and sintered at elevated temperatures, in three zones, up to 1,200°C with programmable controls. Fig. 1 shows the complete fiber extruder assembled in our laboratory. Fig. 2 focuses on the take up spool at the foot of the stand and beneath the heating furnace.

Status of Extruded Fibers

Polycrystalline mullite monofilaments were prepared by the extrusion method. Extrusion was carried out in our new fiber extruder. The orifice for extrusion was varied from 30 μm to 250 μm . Hydrothermally grown KM mullite powder (KM Mullite 101, Kyoritsu Ceramic Materials Co. LTD., Nagoya, Japan) was used as starting powder. The composition of the powder consisted of 60 mol% Al_2O_3 and 40 mol% SiO_2 and the average particle size was 0.77 μm .

PVA and methyl cellulose were used as a binder for extrusion. 15 wt% of PVA (Polyvinyl Alcohol, Airvol 540S, Celanese Chemicals, USA) solution in DI water was prepared as PVA binder. Also 3 wt% of cellulose (Hydroxypropyl Methyl Cellulose, Sigma-Aldrich, Milwaukee, WI) binder was added. The prepared mullite powder was mixed with the organic binder to create a paste for extrusion. The mixed paste was loaded in the fiber extruder and mullite fiber was extruded. The speed of the piston during extrusion was controlled in order to form a uniform diameter of fiber (Fig. 3). The extruded fibers were dried in air and sintered at various temperatures up to 1600 °C. The microstructure of sintered fiber was characterized by scanning electron microscopy (Zeiss DSM960 SEM).




Fig. 1 Photograph of complete Marksman fiber extruder set up at UIUC



Fig. 2. The take up spool of the Marksman fiber extruder




Figure 3 Mullite fiber drawn with the Marksman extruder and sintered at (1600 °C/1hr). This fiber will now be passed through the quadrupole furnace at an optimized temperature and traverse rate.